

# Supermouth®

ultim8™  
SmartBrush System




[supermouth.com](http://supermouth.com)




The Ultim8 Product



**Revolutionary Bristles**  
Patented bristles conform to teeth contours for superior cleaning and reduced sensitivity.



**UV Sanitizing**  
UV sanitizer eliminates 99.9% of germs for a clean brush every time.



**Waterproof**  
Fully waterproof design for use in the shower and easy cleaning.



**Dual Sonic Modes**  
Two sonic modes for customizable cleaning intensity and effectiveness




**Wireless Charging**  
Convenient wireless charging with SmartMag technology for hassle-free power.



**Smart Timer**  
Built-in timer ensures you brush for the dentist-recommended two minutes.



**Auto Drying**  
Integrated dryer keeps your brush dry and hygienic between uses.



**Pressure Sensor**  
Pressure sensor indicator light prevents overbrushing and protects your gums.



The Ultim8 Business

Roadmap to Maximize Passive Revenue with ULTIM8 SmartBrush

ULTIM8 SmartBrush  
Retail \$99.99



ULTIM8 SmartBrush + SmartHub  
Retail \$169.99



ULTIM8 SmartBrush System  
Retail \$199.99



Potential Earning Example:

5000 active patients • 40% retail profit • 10% first-purchase discount for patients

Total active patients 5,000	ULTIM8 SmartBrush Retail \$99.99	ULTIM8 SmartBrush + SmartHub Retail \$169.99	ULTIM8 SmartBrush System Retail \$199.99
10% patients purchase	\$36 × 500 \$18,000	\$61 × 500 \$30,500	\$72 × 500 \$36,000
15% patients purchase	\$36 × 750 \$27,000	\$61 × 750 \$45,750	\$72 × 750 \$54,000
25% patients purchase	\$36 × 1,250 \$45,000	\$61 × 1,250 \$76,250	\$72 × 1,250 \$90,000
50% patients purchase	\$36 × 2,500 \$90,000	\$61 × 2,500 \$152,500	\$72 × 2,500 \$180,000
75% patients purchase	\$36 × 3,750 \$135,000	\$61 × 3,750 \$228,750	\$72 × 3,750 \$270,000

Your Practice's Potential Earnings:

40% retail profit • 10% first-purchase discount for patients  
Fill in the blank with your actual numbers.

Total active patients patient number	ULTIM8 SmartBrush Retail \$99.99	ULTIM8 SmartBrush + SmartHub Retail \$169.99	ULTIM8 SmartBrush System Retail \$199.99
10% patients purchase	\$36 × patient number retail profit	\$61 × patient number retail profit	\$72 × patient number retail profit
15% patients purchase	\$36 × patient number retail profit	\$61 × patient number retail profit	\$72 × patient number retail profit
25% patients purchase	\$36 × patient number retail profit	\$61 × patient number retail profit	\$72 × patient number retail profit
50% patients purchase	\$36 × patient number retail profit	\$61 × patient number retail profit	\$72 × patient number retail profit
75% patients purchase	\$36 × patient number retail profit	\$61 × patient number retail profit	\$72 × patient number retail profit

Earning levels for SuperMouth Pro participants that appear in this publication are examples and should not be construed as typical or average. Income level achievements are dependent upon the individual and practice business skills, personal ambition, time, commitment, activity and demographic factors.



# The Ultim8 Study



UP  
TO

# 40x-60x

MORE EFFECTIVE

## CONCLUSIONS:

The SuperMouth Ultim8 SmartBrush with SuperBristles Pro removes up to 40x more plaque than the flat-headed manual toothbrush. The SuperBristles Ortho cleans up to 60x better around brackets and wires.

## STUDY:

In Vitro Test to Compare the Ability of Toothbrushes to Clean the Interproximal Regions of a Tooth Model.<sup>1-2</sup>

## AIMS:

- Use acellular human teeth to create a tooth model with a realistically sized interproximal<sup>3</sup> space between the teeth. Interproximal space is the hard-to-clean spaces between teeth where cavities tend to occur.
- Use the tooth model to compare the ability of 2 models of toothbrushes to clean the interproximal<sup>3</sup> regions of a tooth model.

## TEST PRODUCTS:

- SuperMouth Ultim8 Sonic SmartBrush with SuperBristles
- Flathead Manual Toothbrush

## CLEANING COMPARISON

Independent Laboratory Testing Report




Figure 1: Example of Interdental Tooth Model.

Study Reference: 229-LAB-IVC-23-001

1. In Vitro Test to Compare the Ability of Manual Toothbrushes to Clean the Interproximal Regions of a Tooth Model. T. Badrock, G. Thomas & G. Withers. Intertek CRS Elm House, Hooton, UK (2024).

2. In vitro cleaning efficacy and resistance to insertion test of interdental brushes, Diana Wolff, Daniela Joerss, Paul Rau, C.E. Dorfer, Clin Oral Invest (2006), 10:297-3

3. Interdental cleaning helps remove debris and interproximal dental plaque, the plaque that collects between two teeth to reduce the likelihood of gum disease and tooth decay. Interproximal area is the area where the tooth meet.

## Learn more



[supermouth.com](https://supermouth.com)